

Angioway™ Manifolds



Product Features



Visible handle makes the flow control easy and accurate.

Ports: one to three

Pressure: 250psi or 500psi

Extension Tube: optional

ON valve and OFF valve are available

Sterilization: ETO

Ordering Information



Catalogue No	Pressure	Side Ports	On/Off Type	Extension Tube
LPJD251N	150psi	stopcock	On	Without
LPJD251F	150psi	stopcock	Off	Without
LPJD252N	250psi	2	On	Without
LPJD252F	250psi	2	Off	Without
LPJD253N	250psi	3	On	Without
LPJD253F	250psi	3	Off	Without
LPJD253NX	250psi	3	On	With
LPJD253FX	250psi	3	Off	With
LPJD503N	500psi	3	On	Without
LPJD503F	500psi	3	Off	Without

Packaged 30 units per box

CE Certificated

Contact Us



If you have any query about the product information, please contact our International Sales & Marketing Department:
international@lepumedical.com

SOURCE: Lepu Medical International Sales & Marketing Department