

Rayband™

Radial Artery Compression Tourniquet



Product Features



Design


- Balloon aligned with puncture site minimizes compression to surrounding tissue
- Adjustable band length provides precise fit for patient comfort
- 15ml or 20ml inflator syringe with Luer taper differentiated with standard syringe

Material

- Transparency provides good visibility for the whole hemostasis procedure
 - Medical polymer material with good biocompatibility and optimal stiffness make no irritation
-

Direction for Use





Tips:

- Confirm no damage of the product before use
- Confirm to inflate balloon only through check valve

Ordering Information



Catalogue No.	Balloon Structure	Bandage Length (cm)	Max Injection Air Volume (ml)	Connector
LPYP15	1	24.5	15	Push-pull
LPYP15L	1	27	15	Push-pull
LPYP15XL	1	29	15	Push-pull
LPYP20	2	24.5	20	Push-pull
LPYP20L	2	27	20	Push-pull
LPYP20XL	2	29	20	Push-pull

YPR20	2	24.5	20	Luer
YPR20L	2	27	20	Luer
YPR20XL	2	29	20	Luer

CE Certificated / FDA Approved

Contact Us



If you have any query about the product information, please contact our International Sales & Marketing Department:
international@lepumedical.com

SOURCE: Lepu Medical International Sales & Marketing Department